# The New Pachette 4

from DGH Technology



Globally Recognized as the

Leading Portable, Desktop Pachymeter

Now Featuring

Exclusive

Automatic Report
Generation Using
Wireless Bluetooth®
Technology













The DGH Pachette Series revolutionized the measurement of corneal thickness with a portable, desktop pachymeter that produces accurate, reliable, and easy to obtain results. The New Pachette 4\* also features an exclusive memory system that stores multiple patient records for automatic report generation using Bluetooth® wireless technology.

\*With Bluetooth® Module installed



# **EASE OF USE**

Turn on power and the Pachette 4 is ready to take measurements. Contact the probe tip to any location on the cornea, and the Pachette 4 will obtain and store up to 25 measurements in a fraction of a second. Audible feedback indicates when a valid measurement is complete. After bilateral measurements have been completed, the measurements can be reviewed on the 16 character by 2 line LCD which provides easy visibility of the current measurement, the average, and the standard deviation of all measurements taken.

## **ACCURACY & REPRODUCIBILITY**

Proven measurement algorithm yields accurate, reproducible measurements in a fraction of a second. The measurement capability of the Pachette 4 is an enhancement of the proven technology that has set the standard for measuring corneal thickness since 1982.





# PERSONALIZED CONFIGURATION MEMORY

When you receive your Pachette 4, the unit is ready to take corneal measurements at initial power up. It is not necessary to set up or configure any parameters. However, the Pachette 4 has been designed to allow the user to personalize the operational mode, number of measurements taken, adjustment of time delays, etc. Once modified, these parameters are permanently stored in its memory.

## **DURABLE PORTABILITY**

The Pachette 4 is designed for performance, portability and convenience, while using limited counter space. The Pachette 4 is encased in a durable rubber holster with a built-in tilt stand. The probe is detachable for easy cleaning and replacement. The probe can be stored in the custom molded probe holder cavity and cord wrap, protecting the probe during transportation or storage. The Pachette 4 has a minimum lifetime of 10 years and an industry leading 3 year warranty.



# **ADDITIONAL FEATURES**

Mapping Mode for obtaining and storing up to 33 actual and biased corneal mapped measurements. These measurements can be obtained at any point on the cornea, including the periphery.

"Off the shelf" batteries: The Pachette 4 is supplied with 2 sets of NiMh rechargeable batteries (4 total AA batteries), or alkaline batteries can also be used. A battery charger is included with each unit.

# Multiple Patient Records Using Bluetooth® Wireless Technology and DGH Connect Software











- Using the DGH Connect Software (PC), enter patient ID information for up to 20 patients and wirelessly transmit to the Pachette 4.
- The Pachette 4 stores all patient ID information in its memory for reference when measuring patients.
- After the measurement sequence is complete, wirelessly transmit all results and patient information directly to a Bluetooth® enabled printer.
  - OR -
- After the measurement sequence is complete, wirelessly transmit all results and patient information back to the PC for report generation and printing.

DGH Technology, Inc. Pachette 4 Ultrasonic Pachymeter SW Ver: 3.0.0

#### PO1 PATIENT INFO

L Name: Smith
F Name: Abraham
ID: 10231
DOB: Jan-01-1971
Gender: Male

#### MEASUREMENT INFO

Date of Meas: Apr-14-2014 Time of Meas: 10:00:36 AM Meas Mode: Cont. Averaging Corneal Vel: 1640 m/s

#### OD MEASUREMENTS

NO.	MEAS(μm)	NO.	MEAS(μm)
1	510	14	510
2	509	15	510
3	509	16	511
4	509	17	510
5	508	18	510
6	509	19	510
7	509	20	510
8	508	21	509
9	509	22	510
10	509	23	509
11	509	24	510
12	509	25	510
13	510		

Average of 25 Meas =  $509 \mu m$ Standard Deviation =  $0.7 \mu m$ 

#### OS MEASUREMENTS

NO.	MEAS(µm)	NO.	MEAS(µm)
1	512	14	513
2	511	15	513
3	512	16	514
4	513	17	513
5	513	18	514
6	513	19	512
7	513	20	514
8	512	21	514
9	514	22	514
10	512	23	512
11	513	24	512
12	512	25	514
13	512		

Average of 25 Meas = 513 μm Standard Deviation = 0.9 μm

# NOTES

Figure A: Bluetooth® enabled printer sample printout (Not shown at actual size)

THE ULTRASOUND SPECIALIST
Serving Eye Care Professionals since 1982

Name: Abraham Smith Pachymetry Report ID: 10231 DGH 555B SN125 DOB: Jan-01-1971 v3.0.0 Gender: Male Apr-14-2014 10:00:36 Notes: Corneal Velocity 1640 m/s OD Average (μm): 509 OS Average (μm): 513 Std Dev (µm): 0.7 Std Dev (µm): 0.9 Meas (μm) Meas (μm) 510 509 511 509 512 509 513 508 513 509 513 509 513 508 512 509 514 509 512 509 513 509 512 510 512 510 513

15 513

18

19 512

22

514

513

514

514

512

Figure B: PC Automatic Report Generation

510

511

510

510

510 509

510 509

510 510

18

19 510

22

(Not shown at actual size)

#### **Specifications**

Measurement Range: Standard Range: 200 microns to 1100 microns,

Flap Option Range: 95 microns to 1100 microns

**Accuracy:** ± 5 microns **Resolution:** 1 micron

**Internal Clock Accuracy:** ±0.002%

**Dimensions:** 6.3" (160 mm) Length x 4.1" (104 mm) Width x 2.2" (56 mm) Depth

Weight: 1.4 lbs (635 grams)

Power Requirement: Two (2) AA batteries

#### **DGH Connect Software: System Requirements**

**Processor:** 32-bit or 64-bit, 2 GHz

Memory: 2 GB RAM

Hard Drive: 1 GB minimum, 100 GB recommended

Ports: USB 2.0

**Display:** 1024 x 768 resolution

Peripherals: Mouse (or Touchpad), Keyboard

Bluetooth® Radio: v2.1 or later\*

\* Software is compatible with USB Bluetooth® adapters/dongles for desktop systems.

# **DGH Connect Software: Compatible Operating Systems**

Microsoft Windows XP, Service Pack 3 or higher (32-bit)

Microsoft Windows Vista, Service Pack 2 or higher (32-bit or 64-bit)

Microsoft Windows 7 (32-bit or 64-bit) Microsoft Windows 8/8.1 (32-bit or 64-bit)

110 Summit Dr., Suite B Exton, PA 19341 Tel. (800) 722-3883 • Fax (610) 594-0390 Email: info@dghkoi.com www.dghkoi.com